



# SKR054-ET

## EMERGENCY TROLLEY

Register code: SKM-A



**SAIKANG**/Medical

Healthier  
your life

# SKR054-ET

## FEATURES AND HIGHLIGHTS


### Technical parameters

- |                         |                |
|-------------------------|----------------|
| ▪ External size (LxWxH) | 750x520x960 mm |
| ▪ Castors               | 100 mm         |



#### One-piece countertop

The one-piece countertop provides more comprehensive protection; the fully guarded handle makes it more comfortable to hold me; the raised grip shape makes the operation of the side-drawing board more comfortable.



#### Decorative gripper

Outer beveled decorative grip, better operating comfort and more three-dimensional shape.



#### Heavy duty construction

Wheels assembled on the same vertical line as four columns. More stable, even the drawer open inside full of medicine/tools.



#### Drawer

Drawer insize big & deeper!



#### Power supply

The power socket is movable and detachable, easy to use, located above the CPR board to convenient the use of oxygen cylinders



#### Inserted oxygen bottle holde

Inserted retractable oxygen bottle holder, belts fixed on trolley main body, space saving



#### Infusion stand

IV pole inserted into main body column.



#### Drawer holder

Card insert from downside, easy keep clean.



#### Defibrillator

The plate surface is made of virgin ABS material. The position of the I,V pole and the defibrillator can be exchanged for easy use.



#### Sorting trash can

Green+yellow sorted trash can.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SAIKANG** / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎ +86 18021231901 · 🕒 24/7  
✉ sales@chinasaikang.com  
🌐 www.saikangmedical.com



Company video



General catalogue